

SINGLE PATIENT USE NIBP CUFFS

- 13 cm

CELatex

4 - 8 cm

CELatex

REV.06 (10.02.2012)

FMT Single Patient Use blood pressure cuffs were designed to reduce the spread of infection through cross-contamination. FMT offers a complete line of single patient use NIBP cuffs from Large Adult size to Neonatal small size. All cuffs are manufactured from a soft, durable, white TPU material which can be quickly and easily wiped The cuff has been manufactured with clean. rounded corners to minimize discomfort.



Designed for single patient use to eliminate cross-contamination

Durable hook-and-loop closure allows repeated repositioning for ease of use

Extra soft and comfortable TPU material is durable and 100% latex-free.

Wide range of sizes to accommodate all patients. Cuffs are available in 12 sizes.

In Neonate sizes color coded printing on each cuff makes it easy to identify cuff size.

Dual or single tube styles to suit the configuration of your blood pressure equipment.

Range guide and index marking on the cuff help to determine the proper size cuff for the patient.

- Adapters available for any monitor on the market
- Cuffs have 20 cm tube/s
- Packaged 10 or 40 per case

FMT #1 to #5 sizes Neonatal cuffs are supplied with a Luer slip connector. Standard male Luer connectors allow FMT neonatal blood pressure cuffs to be adapted to a wide range of most major monitor manufacturers. FMT #6 to #12 sizes cuffs are shipped without connector or with the plastic connectors attached that allow them to be used NIBP with virtually any monitor or sphygmomanometer. FMT #6 to #12 sizes Single Patient Use blood pressure cuffs are available with any of the following connectors in either one- or two-tube versions.

Neonatal

FMT Neocuff

THIS SIDE TO PATIENT

SINGLE PATIENT US

4 7 - 13 cm

CELATE

CELAter



SINGLE PATIENT USE NIBP CUFFS

FMT NeoCuff Single Patient Use Neonatal NIBP Cuffs – Single Tube					
		ff Dimensions L x W	Cuff Range		
REF:	Size	(cm)	(cm)	Color	
FMT-DS/5	#5 Infant	19.1 x 5.6	8.0 - 15.0	Orange 📃	
FMT-DS/4	#4 Neonatal Extra Large	18.1 x 4.8	7.0 - 13.0	Red	
FMT-DS/3	#3 Neonatal Large	14.6 x 4.1	6.0 - 11.0	Green	
FMT-DS/2	#2 Neonatal Medium	11.9 x 3.1	4.0 - 8.0	Light Blue	
FMT-DS/1	#1 Neonatal Small	10.0 x 2.5	3.0 - 6.0	Navy Blue	

FMT NeoCuff Single Patient Use Neonatal NIBP Cuffs – Double Tube				
	Cuf	f Dimensions L x W	Cuff Range	
REF:	Size	(cm)	(cm)	Color
FMT-DD/5	#5 Infant	19.1 x 5.6	8.0 – 15.0	Orange
FMT-DD/4	#4 Neonatal Extra Large	18.1 x 4.8	7.0 - 13.0	Red 📃
FMT-DD/3	#3 Neonatal Large	14.6 x 4.1	6.0 - 11.0	Green
FMT-DD/2	#2 Neonatal Medium	11.9x 3.1	4.0 - 8.0	Light Blue
FMT-DD/1	#1 Neonatal Small	10.0 x 2.5	3.0 - 6.0	Navy Blue

FMT Single Patient Use NIBP Cuffs – Single Tube					
		Cuff Dimensions L x W	Cuff Range		
REF:	Size	(cm)	(cm)	Color	
FMT-DS/12	#12 Large Adult	63.0 x 17.0	33.0 - 47.0	White	
FMT-DS/11	#11 Adult	51.0 × 14.0	25.3 - 34.3	White	
FMT-DS/10	#10 Small Adult	40.0 x 11.4	20.0 - 27.0	White	
FMT-DS/9	#9 Pediatric	38.0 x 8.6	15.8 - 21.3	White	
FMT-DS/8	#8 Small Pediatric	30.0 x 7.2	12.4 - 16.8	White	
FMT-DS/7	#7 Infant	25.0 x 6.0	9.8 - 13.3	White	
FMT-DS/6	#6 Neonate	20.0 x 5.0	7.7 - 10.5	White	

FMT Single Patient Use NIBP Cuffs – Double Tube					
		Cuff Dimensions L x W	Cuff Range		
REF:	Size	(cm)	(cm)	Color	
FMT-DD/12	#12 Large Adult	63.0 x 17.0	33.0 - 47.0	White	
FMT-DD/11	#11 Adult	51.0 x 14.0	25.3 - 34.3	White	
FMT-DD/10	#10 Small Adult	40.0 x 11.4	20.0 - 27.0	White	
FMT-DD/9	#9 Pediatric	38.0 x 8.6	15.8 - 21.3	White	
FMT-DD/8	#8 Small Pediatric	30.0 x 7.2	12.4 - 16.8	White	
FMT-DD/7	#7 Infant	25.0 x 6.0	9.8 - 13.3	White	
FMT-DD/6	#6 Neonate	20.0 x 5.0	7.7 - 10.5	White	



Please add the connector numbers at the end of the cuff reference number for the complete cuff part number. Example: FMT-DS/12BP12, FMT-DD/10BP0202



FMT[®] is registered trademark of Metko Ltd. Other products, logos and company names mentioned herein maybe trademarks of their respective owners and are not to be taken as an endorsement or affiliation with FMT[®]. The information in this instruction insert has been carefully checked and it is believed to be accurate. In the interest of continued product development, METKO reserves the right to make changes and improvements to this insert and the product it described any time, without notice or obligation.



HEADQUARTER & FACTORY -

İvedik O.S.B. Ağaç İşleri Sanayi Sit. 1354. Cad. 1358.Sok. No:9 06378 Yenimahalle

ANKARA / TÜRKİYE

Tel (pbx) : +90 312 387 12 46 Fax : +90 312 387 12 51

BRANCH

Atatürk Mah. Sedef Cad. 38 Ada Ata Plaza 3-5 Blokları Kat:2 No:57 34758 Ataşehir

İSTANBUL / TÜRKİYE

Tel (pbx) : +90 216 455 91 21 Fax : +90 216 455 48 24

INTERNET

http://www.metkoltd.com metko@metkoltd.com